

**Laminar Flow Clean Benches, Horizontal and Vertical** 

The Leading Solution for Research Laboratories

More product variants and sizes to suit more applications and user preferences

### **KEY BENEFITS:**

- Energy-Efficient DC ECM Blower. ONLY 155 Watts!! (for 4 ft size) and with Standby Mode for more energy saving
- World's Leading Low Noise Laminar Flow Cabinet
- Stable and Self-compensating Flow
- Isocide™ Powder-coated
- ULPA filter (up to ISO Class 3)
- With Real-time Airflow Monitoring



### **Airflow Sensor**

- Monitors real-time airflow for safety
- Alerts the user if the airflow is insufficient

ESCO



#### **Rocker Switches and Pressure Gauge**

- **■** Easy-to-use switches
- Displays filter loading status
- Standby mode switch for more energy savings



### Optional Sentinel™ Gold Microprocessor Controller

- Displays all safety information on one screen
- Standby mode switch to further save energy
- Programmable UV timer to extend UV lamp life

### **Ergonomic Control Panel Location**

- Located on the center and angled down
- ADA-Compliant



### Work Top \_\_\_\_

 The spill-retaining work top design with a recessed central area contains accidental liquid spills.

I H

#### **Guide to Models**

	<del></del>							_			_		_		
									_		7				Ţ
1st Placeholo	ler	2nd Placeho	older	3rd Placehol	der	4th Placeh	nolder	5th Placel	nolder	6th Placehold	er	7th Placeho	lder	8th Placehold	ler
Product Lin	ie	Flow		Side Wall		Widt	h	Internal F	leight	Control		Windov	v	Electrical	
Laminar Flow	L	Horizontal	н	Glass	*G	3 feet	3	2.5 feet	В	Simple Switches Standard	S	Fixed	F	115 VAC, 50 / 60 Hz	9
				Stainless Steel	S	4 feet	4	3 feet	***C	Sentinel™ Gold Optional	G				
						5 feet	5								
						6 feet	6								
						**8 feet	8								

<sup>\*</sup> Glass side is only available for 2 ½ feet (B) internal height Example Code: LHG-4BS-F9
Specifications: Laminar flow cabinet, Horizontal Flow, Glass sided walls, 4 feet width, 2.5 feet internal height, Simple Switches controller, Fixed window and 115 VAC, 50/60 Hz

2


<sup>\*\*8</sup> feet width is only available for 3 feet (C) internal height Example Code: LHS-8CS-F9
Specifications: Laminar flow cabinet, Horizontal Flow, Stainless steel sided walls, 8 feet width, 3 feet internal height, Simple Switches controller, Fixed window and 115 VAC, 50/60 Hz


<sup>\*\*\*3</sup> feet internal height is only available for stainless steel sided walls AND 4 feet, 6 feet or 8 feet width Example Code: LHS-4CS-F9
Specifications: Laminar flow cabinet, Horizontal Flow, Stainless steel sided walls, 4 feet width, 3 feet internal height, Simple Switches controller, Fixed window and 115 VAC, 50/60 Hz

## **Energy Efficient DC ECM Blower**

- Powered by latest generation DC ECM, that is more efficient than legacy ECM and VFD motors
- 70% Energy savings compared to AC motor
- Stable airflow, despite building voltage fluctuations & filter loading
- Standby mode to further reduce power consumption by 60%







### **Quiet Operation**

- World's leading low noise laminar flow cabinet
- Comfortable low noise at 55 dBA


\*LHG-4B in open field condition

Reduces fatigue and improves work concentration



#### **ULPA Filter**

- 10x Filtration efficiency of HEPA filter Creates ISO Class 3 work zone instead of industry-standard ISO Class 5
- 10x Cleaner work zone than cabinets with HEPA filter




#### Airstream® Horizontal Laminar Flow Clean Bench, Model LHG-\_B Available in 0.9, 1.2, 1.5, and 1.8 meter models (3', 4', 5', and 6')

111 🔾

#### **ISOCIDE™** Powder Coat

- Silver-ion impregnated powder coat
- Inhibits microbial growth to improve safety



	Cabinet Performance	Air Quality	Filtration	Electrical Safety
Standards Compliance	IEST-RP-CC002.2, Worldwide	ISO 14644.1, Class 4, Worldwide AS 1386 Class 2.5, Australia JIS B9920, Class 4, Japan	EN-1822 (H14), Europe IEST-RP-CC001.3, Worldwide IEST-RP-CC007, Worldwide IEST-RP-CC034.1, Worldwide	IEC61010-1, Worldwide EN 61010-1, Europe UL61010-1, USA CAN/CSA-22.2, No.61010-1



# Optional Airflow Sensor

- Monitors real-time airflow for safety
- Alerts the user if airflow is insufficient

ESCO



90.0

#### **Rocker Switches and Pressure Gauge**

- **■** Easy-to-use switches
- Displays filter loading status
- Standby mode switch for more energy savings



### **Optional Sentinel™ Gold Microprocessor Controller**

- Displays all safety information on one screen
- Standby mode switch to further save energy
- Programmable UV timer to extend UV lamp life

### **Ergonomic Control Panel Location**

- Located on the center and angled down
- ADA-Compliant



### Work Top

 Spill-retaining work top design with a recessed central area contains accidental liquid spills.






	L V										-				
									$\overline{}$		$\overline{\mathcal{L}}$				
1st Placehold  Product Lin		2nd Placeho	older	3rd Placehold Side Wall		4th Placel		5th Placeh		6th Placehold <b>Control</b>	er	7th Placeho <b>Windov</b>		8th Placehold <b>Electrical</b>	
Product Em	e	FIOW		Side Wall		vviac	.11	IIIterriai r	leigiit			VVIIIdov	V		
Laminar Flow	L	Vertical	V	Glass	*G	3 feet	3	2 ¼ feet	A	Simple Switches Standard	S	Fixed	F	115 VAC, 50 / 60 Hz	9
				Stainless Steel	S	4 feet	4	3 feet	***C	Sentinel <sup>™</sup> Gold Optional	G				
						5 feet	5								
						6 feet	6								
						**8 feet	8								


<sup>\*</sup> Glass side is only available for 2 ¼ feet (A) internal height Example Code: LVG-4AS-F9
Specifications: Laminar flow cabinet, Vertical Flow, Glass sided walls, 4 feet width, 2 feet internal height, Simple Switches controller, Fixed window and 115 VAC, 50/60 Hz

### **Energy Efficient DC ECM Blower**

- Powered by latest generation DC ECM, that is more efficient than legacy ECM and VFD motors
- 70% Energy savings compared to AC motor
- Stable airflow, despite building voltage fluctuations & filter loading
- Standby mode to further reduce power consumption by 60%







### **Quiet Operation**

- World's leading low noise laminar flow cabinet
- Comfortable low noise at 52 dBA


\*LVG-4A in open field condition

Reduces fatigue and improves work concentration



#### **ULPA Filter**

- 10x Filtration efficiency of HEPA filter Creates ISO Class 3 work zone instead of industry-standard ISO Class 5
- 10x Cleaner work zone than cabinets with HEPA filter



#### **ISOCIDE™** Powder Coat

- Silver-ion impregnated powder coat
- Inhibits microbial growth to improve safety



Airstream® Vertical Laminar Flow Clean Bench Model LVG-\_A Available in 0.9, 1.2, 1.5 , and 1.8 meter models (3', 4', 5', and 6').

10


	Cabinet Performance	Air Quality	Filtration	Electrical Safety
Standards Compliance	IEST-RP-CC002.2, Worldwide	ISO 14644.1, Class 4, Worldwide AS 1386 Class 2.5, Australia JIS B9920, Class 4, Japan	EN-1822 (H14), Europe IEST-RP-CC001.3, Worldwide IEST-RP-CC007, Worldwide IEST-RP-CC034.1, Worldwide	IEC61010-1, Worldwide EN 61010-1, Europe UL61010-1, USA CAN/CSA-22.2, No.61010-1




### **Cabinet Filtration System**


- Room air is taken in from the top of the clean bench through a disposable pre-filter with 85% arrestance; this serves to trap larger particles and increase the life of the main filter.
- Air is forced evenly across the ULPA/H14 filter(s); the result is a stream of clean laminar air within the workzone of the clean bench; this dilutes and flushes all airborne contaminants from the interior.
- A nominal filter face velocity of 0.45 m/s or 90 fpm ensures that there is a sufficient number of air changes within the enclosed area of the clean bench to maintain cleanliness.
- The purified air travels across the internal work zone of the clean bench in a horizontal, unidirectional stream and leaves the main work chamber across the entire open front of the clean bench.

### **Airflow Diagram**



- ULPA / H14 filtered air
- Room air / Inflow air





- Microprocessor Controller / Simple Switches Control System
- 2. UV Light Retrofit Kit Provision
- 4. Spill-retaining Stainless Steel Work Top
- 5 Pre-filter
- 6. DC ECM Blower

- 7. Fluorescent Lamp
- 8. Electrical Outlet Provision (maximum of 2)

5

6

8

10

11

12

13

- 9. IV Bar Retrofit Kit Provision
- 10. Air Flow Sensor
- LH\_-A: 1118 mm (44.0") / LH\_-B: 1270 mm (50.0") / LH\_-C: 1422 (56.0" mm) (22.5") (28.5") (34.5") 573 725 877 神神で **H**H 631 mm (24.8") for LHG Models 620 mm (24.4") for LHS Models LH\_-A: 795 mm (31.3") LH\_-B: 788 mm (31.0") LH\_-C: 782 mm (30.8")
  - 11. Service Fixture Retrofit Kit Provision (2 on each side)
  - 12. Tempered Glass Side Walls (for LHG Variant) Stainless Steel Side Walls (for LHS Variant)
  - 13. ULPA / H14 Filter

	Esco Horizontal and Vertical Laminar Flow Clean Benches Accessories									
Model	3 feet Width	4 feet Width	4 feet Width 5 feet Width 6 feet Width							
UV Lamp	UV-15A-L 5170251		UV-30A-L 5170255		UV-30A-L *2					
Electrical Outlet			EO-GFCI 5170071							
IV Bar Horizontal Flow	IV-890 5170608	IV-1195 5170609	IV-1500 5170610	IV-1805 5170611	IV-2375 5170662					
IV Bar Vertical Flow	IV-960 5170603	IV-1265 5170604 IV-1570 5170605 IV-1875 5170606 IV-2445 51								
PVC Arm Rest	ME-W-REST 5170127									
Height-adjustable Lab Chair (Blue)	ME-LD-AR360 1150006									
Ergonomic Foot Rest			FT-REST 5170073							
Support Stand Telescoping with Castors (28" to 34")	STC-3A0 5130055	STC-4A0 5130056	STC-5A0 5130057	STC-6A0 5130058	STC-8A0 5130059					
Support Stand Telescoping with Leveling Feet (28" to 34")	STL-3A0 5130050	STL-4A0 5130051	STL-5A0 5130052	STL-6A0 5130053	STL-8A0 5130054					
Support Stand with Leveling Feet (28")	SAL-3A0 Gen2 5130170	SAL-4A0 Gen2 SAL-5A0 Gen2 5130134 5130171		SAL-6A0 Gen2 5130172	SAL-8A0 Gen2 5131124					
Support Stand with Leveling Feet (34")	SAL-3B0 Gen2 5130174	SAL-4B0 Gen2 5130175	SAL-5B0 Gen2 5130176	SAL-6B0 Gen2 5130177	SAL-8B0 Gen2 5131125					


6




### **Cabinet Filtration System**


- Room air is taken in from the top of the clean bench through a disposable pre-filter with 85% arrestance; this serves to trap larger particles and increase the life of the main filter.
- Air is forced evenly across the ULPA/H14 filter(s); the result is a stream of clean laminar air within the workzone of the clean bench; this dilutes and flushes all airborne contaminants from the interior
- A nominal filter face velocity of 0.45 m/s or 90 fpm ensures that there is a sufficient number of air changes within the enclosed area of the clean bench to maintain cleanliness.
- The purified air travels across the working zone of the clean bench in a vertical, unidirectional stream and leaves the main work chamber across the entire open front of the clean bench and through Auto-Purge™ slots at the back wall of the work zone which are designed to eliminate air turbulence and the possibility of dead air corners in the work zone.

### **Airflow Diagram**



- Room air / Inflow air





- 1. Esco Sentinel™ Gold Microprocessor Control System / Simple Switches Control System
- 2. Fixed Sash (LV\_-F variant) (Optional) Manual Sliding Sash (LV\_-S variant)
- 3. Optional Front Cover
- 4. Stainless Steel Work Top
- 5 Pre-filter

- 6. DC ECM Blower (Self-compensating and lower noise)
- 7. ULPA filter
- 8. IV Bar Retrofit Kit Provision
- 9. Fluorescent Lamp
- 10. Tempered Glass Side Walls (for LVG Variant) Stainless Steel Side Walls (for LVS Variant)
- 11. Airflow Sensor
- 12. UV Lamp
- 13. Service Fixture Retrofit Kit Provision (2 holes on each side)
- 14. Electrical Outlet Retrofit Kit provison (1 outlet on each side)

#### Esco Horizontal and Vertical Laminar Flow Clean Benches Front Cover (For UV Interlock System) 3 feet 4 feet 5 feet 8 feet 6 feet FC-LHG / LHS-3A FC-LHG / LHS-4A FC-LHG / LHS-5A FC-LHG / LHS-6A 5170601 5170602 5170585 5170586 FC-LHG / LHS-3B FC-LHG / LHS-4B FC-LHG / LHS-5B FC-LHG / LHS-6B LHG / LHS 5170589 5170590 5170587 5170588 FC-LHG / LHS-3C FC-LHG / LHS-4C FC-LHG / LHS-6C FC-LHG / LHS-8C 5170627 5170591 5170592 5170593 FC-LVG / LVS-3A FC-LVG / LVS-4A FC-LVG / LVS-5A FC-LVG / LVS-6A 5170595 5170584 5170596 5170597 LVG / LVS\* FC-LVG / LVS-4C FC-LVG / LVS-5C FC-LVG / LVS-6C FC-LVG / LVS-8C 5170600 5170598 5170061 5170599

<sup>\*</sup> For fixed sash models only (LVS- F and LVG- F)

### General Specifications, Airstream® Horizontal Laminar Flow Clean Benches, B-Height (Interior Height: 2 ½ ft / 0.8 m)

Glass Side, Simple Switches		LHG-3BS-F9 2120705	LHG-4BS-F9 2120716	LHG-5BS-F9 2120717	LHG-6BS-F9 2120718				
Stainless Steel Side, Simple Switches		LHS-3BS-F9 2120661	LHS-4BS-F9 2120663	LHS-5BS-F9 2120665	LHS-6BS-F9 2120667				
Stainless Steel Side,	Microprocessor	LHS-3BG-F9 2120503	LHS-4BG-F9 2120505	LHS-5BG-F9 2120507	LHS-6BG-F9 2120509				
Nominal Size		0.9 meter (3')	1.2 meter (4')	1.5 meter (5')	1.8 meter (6')				
External Dimensions (W x D x H)	Without Base Stand	1035 x 788 x 1270 mm (40.8" x 31.0" x 50.0")	1340 x 788 x 1270 mm (52.8" x 31.0" x 50.0")	1645 x 788 x 1270 mm (64.8" x 31.0" x 50.0")	1950 x 788 x 1270 mm (76.8" x 31.0" x 50.0")				
Internal Work Area, (W x D x H)*	. Dimensions	898 x 620 x 725 mm (35.4" x 24.4" x 28.5")	1203 x 620 x 725 mm (47.4" x 24.4" x 28.5")	1508 x 620 x 725 mm (59.4" x 24.4" x 28.5")	1813 x 620 x 725 mm (71.4" x 24.4" x 28.5")				
Internal Work Area,	Space	0.5 m² (5.4 ft²)	0.7 m² (7.5 ft²)	0.9 m² (9.6 ft²)	1.0 m <sup>2</sup> (10.8 ft <sup>2</sup> )				
Average Airflow Ve	locity		0.45 m/s (90 fpm	a) at initial setpoint					
Air Volume		1055 m³/hr (621 cfm)	1413 m³/hr (832 cfm)	1771 m³/hr (1042 cfm)	2129 m³/hr (1253 cfm)				
ULPA Filter Typical E	Efficiency	> 99.999% at particle size between 0.1 to 0.3µm							
Sound Emission per	IEST-RP-CC002.2*	52.8 dBA	55.4 dBA	58.0 dBA	59.6 dBA				
Fluorescent Lamp Ir at Zero Ambient	ntensity	1279 Lux (119 foot candles)	1394 Lux (130 foot candles)	1204 Lux (112 foot candles)	1494 Lux (139 foot candles)				
	Main Body	1.2 mm (0.05") 18-gauge electro-galvanized steel with white oven-baked epoxy-polyster powder-coated finish.							
Cabinet	Work Zone	1.2 mm (0.05") 18-gauge stainless steel, grade 304, with 4B finish.							
Construction	LHG Side Walls	UV-absorbing tempered glass, 5 mm (0.2"), colorless and transparent							
	LHS Side Walls	Interior: 1.2 mm (0.05") 18-gauge stainless steel, grade 304, with 4B finish.  Exterior: 1.2 mm (0.05") 18 gauge electro-galvanized steel with white oven-baked epoxy-polyester powder-coated finish.							
	Cabinet Full Load Amps (FLA)	12.3 A	12.4 A	12.4 A	12.5 A and 6.5 A (2 Power Inlet)				
Electrical Power Rating 9:	Optional Outlets (FLA)	6 A							
110-130 VAC, 50 / 60Hz, 1 phase**	Cabinet Nominal Power (W)	154	188	273	343				
	Heat Rejected, BTU per Hour	525	641	932	1170				
Net Weight***		108 Kg (238 lbs)	156 Kg (344 lbs)	180 Kg (397 lbs)	228 Kg (503 lbs)				
Shipping Weight***		140 Kg (308 lbs)	211 Kg (465 lbs)	237 Kg (522 lbs)	293 Kg (646 lbs)				
Shipping Dimensions, Maximum (W x D x H)***		1120 x 900 x 1590 mm (44" x 35" x 62")	1400 x 900 x 1590 mm (55" x 35" x 62")	1720 x 900 x 1590 mm (68" x 35" x 62")	2200 x 900 x 1590 mm (87" x 35" x 62")				
Shipping Volume, M	laximum***	1.6 m³ (56.6 ft³)	2.0 m³ (70.6 ft³)	2.5 m³ (87 ft³)	3.3 m³ (118.6 ft³)				

<sup>\*</sup>Noise reading in open field condition / anechoic chamber. Noise reading in normal room varies by room size, layout, and background noise, but may reach roughly 3 dBA above these values.

\*\*Additional voltages may be available; contact Esco for ordering information.

\*\*\*Clean bench only; excludes optional stand.

### General Specifications, Airstream® Horizontal Laminar Flow Clean Benches, C-Height (Interior Height: 3 ft / 0.9 m)

Stainless Steel Side, Simple Switches		LHS-4CS-F9 2120675	LHS-6CS-F9 2120679	LHS-8CS-F9 2120681		
Nominal Size		1.2 meter (4')	1.8 meter (6')	2.4 meter (8')		
External Dimensions (W x D x H) Without Base Stand		1340 x 788 x 1422 mm (52.8" x 31.0" x 56.0")	1950 x 788 x 1422 mm (76.8" x 31.0" x 56.0")	2520 x 788 x 1422 mm (99.2" x 31.0" x 56.0")		
nternal Work Area, (W x D x H)*	Dimensions	1203 x 620 x 877 mm (47.4" x 24.4" x 34.5")	1813 x 620 x 877 mm (71.4" x 24.4" x 34.5")	2383 x 620 x 877 mm (93.8" x 24.4" x 34.5")		
nternal Work Area,	Space	0.7 m <sup>2</sup> (7.5 ft <sup>2</sup> )	1.0 m <sup>2</sup> (10.8 ft <sup>2</sup> )	1.4 m <sup>2</sup> (15.0 ft <sup>2</sup> )		
Average Airflow Ve	elocity		0.45 m/s (90 fpm) at initial set point			
Air Volume		1709 m³/hr (1006 cfm)	2576 m³/hr (1516 cfm)	3385 m³/hr (1992 cfm)		
JLPA Filter Typical I	Efficiency	>	· 99.999% at particle size between 0.1 to 0.3μn	n		
ound Emission per	r IEST-RP-CC002.2*	56.4 dBA	59.4 dBA	62.3 dBA		
Fluorescent Lamp Ir Ambient	ntensity At Zero	1304 Lux (121 foot candles)	1136 Lux (106 foot candles)			
	Main Body	1.2 mm (0.05") 18-gauge electr	ro-galvanized steel with white oven-baked epoxy-polyster powder-coated finish.			
Cabinet Construction	Work Zone	1.2 mm (0	0.05") 18-gauge stainless steel, grade 304, with	4B finish.		
	LHS Side Walls			gauge stainless steel, grade 304, with 4B finish. ed steel with white oven-baked epoxy-polyester powder-coated finish.		
	Cabinet Full Load Amps (FLA)	12.3 A	12.5 A and 6.5 A (2 Power Inlet)	13 A and 7 A (2 Power Inlet)		
lectrical lower Rating 9:	Optional Outlets (FLA)		6 A			
10-130 VAC, 50 / 60Hz, phase**	Cabinet Nominal Power (W)	217	400	432		
	Heat Rejected, BTU per Hour	740	1365	1481		
let Weight***		167 Kg (368 lbs)	236 Kg (520 lbs)	317 Kg (699 lbs)		
Shipping Weight***		222 Kg (489 lbs)	301 Kg (663 lbs)	397 Kg (875 lbs)		
hipping Dimension W x D x H)***	ns, Maximum	1400 x 900 x 1590 mm (55" x 35" x 62")	2200 x 900 x 1590 mm (87" x 35" x 62")	2720 x 950 x 1590 mm (107" x 37" x 62")		
hipping Volume, N	Maximum***	2.0 m³ (70.6 ft³)	3.3 m³ (118.6 ft³)	4.1 m³ (144.8 ft³)		



<sup>\*</sup>Noise reading in open field condition / anechoic chamber. Noise reading in normal room varies by room size, layout, and background noise, but may reach roughly 3 dBA above these values.

\*\*Additional voltages may be available; contact Esco for ordering information.

\*\*\*Clean bench only; excludes optional stand.

### General Specifications, Airstream® Vertical Laminar Flow Clean Benches, A-Height (Interior Height: 2 ¼ ft / 0.7 m)

Glass Side, Microprocessor		LVG-3AG-F9 2120435	LVG-4AG-F9 2120437	LVG-5AG-F9 2120439	LVG-6AG-F9 2120441				
Stainless Steel Side,	Microprocessor	LVS-3AG-F9 2120443	LVS-4AG-F9 2120445	LVS-5AG-F9 2120447	LVS-6AG-F9 2120449				
Nominal Size		0.9 meter (3')	1.2 meter (4')	1.5 meter (5')	1.8 meter (6')				
External Dimensions (W x D x H)	Without Base Stand	1035 x 824 x 1270 mm (40.8" x 32.4" x 50.0")	1340 x 824 x 1270 mm (52.8" x 32.4" x 50.0")	1645 x 824 x 1270 mm (64.8" x 32.4" x 50.0")	1950 x 824 x 1270 mm (76.8" x 32.4" x 50.0")				
Internal Work Area, Dimensions	LVG Models	965 x 739 x 689 mm (38.0" x 29.1" x 27.1")	1270 x 739 x 689 mm (50.0" x 29.1" x 27.1")	1575 x 739 x 689 mm (62.0" x 29.1" x 27.1")	1880 x 739 x 689 mm (74.0" x 29.1" x 27.1")				
(W x D x H)*	LVS Models	965 x 739 x 678 mm (38.0" x 29.1" x 26.7")	1270 x 739 x 678 mm 1575 x 739 x 678 mr ) (50.0" x 29.1" x 26.7") (62.0" x 29.1" x 26.7		1880 x 739 x 678 mm (74.0" x 29.1" x 26.7")				
Internal Work Area, S	Брасе	0.6 m² (6.5 ft²)	0.8 m² (8.6 ft²)	1.0 m² (10.7 ft²)	1.3 m² (14.0 ft²)				
Average Airflow Vel	ocity		0.45 m/s (90 fpm	) at initial setpoint					
Air Volume		1117 m³/hr (657 cfm)	1471 m³/hr (866 cfm)	1824 m³/hr (1074 cfm)	2177 m³/hr (1281 cfm)				
ULPA Filter Typical E	fficiency		> 99.999% at particle siz	re between 0.1 to 0.3µm					
Sound Emission per	IEST-RP-CC002.2*	51.6 dBA	52.4 dBA	55.6 dBA	57.6 dBA				
Fluorescent Lamp In at Zero Ambient	tensity	980 Lux (94 foot candles)	904 Lux (84 foot candles)	894 Lux (83 foot candles)	1062 Lux (99 foot candles)				
	Main Body	1.2 mm (0.05") 18-gauge electro-galvanized steel with white oven-baked epoxy-polyster powder-coated finish.							
Cabinet	Work Zone	1.2 mm (0.05") 18-gauge stainless steel, grade 304, with 4B finish.							
Construction	LVG Side Walls	UV-absorbing tempered glass, 5 mm (0.2"), colorless and transparent							
	LVS Side Walls	Interior: 1.2 mm (0.05") 18-gauge stainless steel, grade 304, with 4B finish.  Exterior: 1.2 mm (0.05") 18 gauge electro-galvanized steel with white oven-baked epoxy-polyester powder-coated finish.							
	Cabinet Full Load Amps (FLA)	13 A	13.5 A	14.3 A	14.5 A				
Electrical Power Rating 9: 110-130 VAC,	Optional Outlets (FLA)	6 A							
50 / 60Hz, 1 phase**	Cabinet Nominal Power (W)	132	155	204	264				
	Heat Rejected, BTU per Hour	450	529	696	908				
Net Weight***		135 Kg (298 lbs)	158 Kg (348 lbs)	199 Kg (438 lbs)	208 Kg (459 lbs)				
Shipping Weight***		167 Kg (368 lbs)	202 Kg (445 lbs)	256 Kg (564 lbs)	273 Kg (602 lbs)				
Shipping Dimension (W x D x H)***	s, Maximum	1120 x 900 x 1590 mm (44" x 35" x 62")	1400 x 900 x 1590 mm (55" x 35" x 62")	1720 x 900 x 1590 mm (68" x 35" x 62")	2200 x 900 x 1590 mm (87" x 35" x 62")				
Shipping Volume, M	aximum***	1.6 m³ (56.6 ft³)	2.0 m³ (70.6 ft³)	2.5 m³ (88.2 ft³)	3.3 m³ (116.5 ft³)				

<sup>\*</sup>Noise reading in open field condition / anechoic chamber. Noise reading in normal room varies by room size, layout, and background noise, but may reach roughly 3 dBA above these values.

<sup>\*\*</sup>Additional voltages may be available; contact Esco for ordering information.
\*\*\*Clean bench only; excludes optional stand.

### General Specifications, Airstream® Vertical Laminar Flow Clean Benches, C-Height (Interior Height: 3 ft / 0.9 m)

Stainless Side, Simp	e Switches	LVS-4CS-F9 2120710	LVS-5CS-F9 2120720	LVS-6CS-F9 2120711	LVS-8CS-F9 2120689				
Nominal Size		1.2 meter (4')	1.5 meter (5')	1.8 meter (6')	2.4 meters (8')				
External Dimensions (W x D x H)	Without Base Stand	1340 x 814 x 1645 mm (52.8" x 32.0" x 64.8")	1645 x 814 x 1645 mm (64.8" x 32.0" x 64.8")	1950 x 814 x 1645 mm (76.8" x 32.0" x 64.8")	2520 x 814 x 1645 mm (99.2" x 32.0" x 64.8")				
Internal Work Area, (W x D x H)*	Dimensions	1270 x 739 x 905 mm (50.0" x 29.1" x 35.6")	1575 x 739 x 905 (62.0" x 29.1" x 35.6")	1880 x 739 x 905 mm (74.0" x 29.1" x 35.6")	2450 x 739 x 905 mm (96.5" x 29.1" x 35.6")				
Internal Work Area,	Space	0.8 m² (8.6 ft²)	1.0 m² (10.7 ft²)	1.3 m <sup>2</sup> (14.0 ft <sup>2</sup> )	1.6 m² (17.2 ft²)				
Average Airflow Ve	locity		0.45 m/s (90 fpm	) at initial setpoint					
Air Volume		1432 m³/hr (843 cfm)	1776 m³/hr (1045 cfm)	2120 m³/hr (1248 cfm)	2762 m³/hr (1626 cfm)				
ULPA Filter Typical E	Efficiency		> 99.999% at particle siz	te between 0.1 to 0.3µm					
Sound Emission per	IEST-RP-CC002.2*	52.6 dBA	53.0 dBA	53.2 dBA	59.8 dBA				
Fluorescent Lamp Intensity At Zero Ambient		904 Lux (84 foot candles)	894 Lux (83 foot candles)	1062 Lux (99 foot candles)	1100 Lux (8102 foot candles)				
	Main Body	1.2 mm (0.05") 18-gauge electro-galvanized steel with white oven-baked epoxy-polyster powder-coated finish.							
Cabinet Construction	Work Zone	1.2 mm (0.05") 18-gauge stainless steel, grade 304, with 4B finish.							
	LVS Side Walls		Interior: 1.2 mm (0.05") 18-gauge stainless steel, grade 304, with 4B finish.  1.2 mm (0.05") 18 gauge electro-galvanized steel with white oven-baked epoxy-polyester powder-coated finish.						
	Cabinet Full Load Amps (FLA)	13.5 A	14.3 A	14.5 A and 8.5 A (2 Power Inlet)	15 A and 9 A (2 Power Inlet)				
Electrical Power Rating 9:	Optional Outlets (FLA)		5	A					
110-130 VAC, 50 / 60Hz, 1 phase**	Cabinet Nominal Power (W)	217	304	400	434				
	Heat Rejected, BTU per Hour	740	1037	1365	1481				
Net Weight***		194 Kg (428 lbs)	244 Kg (538 lbs)	255 Kg (560 lbs)	352 Kg (776 lbs)				
Shipping Weight***		243 Kg (536 lbs)	308 Kg (679 lbs)	330 Kg (661 lbs)	455 Kg (1003 lbs)				
Shipping Dimensions, Maximum (W x D x H)***		1490 x 900 x1750 mm (59" x 35" x 69")	1790 x 900 x 1750 mm (70" x 35" x 69")	2070 x 900 x 1750 mm (82" x 35" x 69")	2720 x 950 x 1750 mm (107" x 37" x 69")				
Shipping Volume, M	laximum***	2.3 m³ (81 ft³)	2.8 m³ (99.6 ft³)	3.3 m³ (116.5 ft³)	4.5 m³ (158.9 ft³)				



<sup>\*</sup>Noise reading in open field condition / anechoic chamber. Noise reading in normal room varies by room size, layout, and background noise, but may reach roughly 3 dBA above these values.

<sup>\*\*</sup>Additional voltages may be available; contact Esco for ordering information.
\*\*\*Clean bench only; excludes optional stand.

### **ESCO LIFESCIENCES GROUP**

42 LOCATIONS IN 21 COUNTRIES ALL OVER THE WORLD



Join us on Social Media and Download our Apps!





















































Esco Micro Pte. Ltd. • 21 Changi South Street 1 • Singapore 486 777 Tel +65 6542 0833 • Fax +65 6542 6920 • mail@escolifesciences.com www.escolifesciences.com

Esco Technologies, Inc. • 903 Sheehy Drive, Suite F, Horsham, PA 19044, USA Tel: +1 215-441-9661 • Fax 484-698-7757 eti.admin@escolifesciences.com

Esco Lifesciences Group Offices: Bangladesh | China | Denmark | Germany | Hong Kong | India | Indonesia | Italy | Japan | Lithuania | Malaysia | Myanmar | Philippines | Russia | Singapore | South Africa | South Korea | Taiwan | Thailand | UAE | UK | USA | Vietnam

