

ESCO
LIFESCIENCES GROUP



Labculture® G4

Class II Type A2 Biological Safety Cabinets

*The Most Advanced Energy-efficient, Safe,
and Ergonomic Biosafety Cabinet in the World*




LABCULTURE® G4 (LA2 G4) CLASS II TYPE A2 BIOSAFETY CABINET



USB Port and Zero Volt Relay Contact

- USB Port to send operational information to Building Management System (BMS)
- Zero Volt Relay Contact to turn ON/OFF exhaust blower and signal the building alarm



Airflow Sensor

- Monitors real-time airflow for safety
- Alert the user if airflow is insufficient



Centurion 7" Capacitive Touchscreen Controller

- Displays all safety information on one large screen
- Shows cabinet parameters with intuitive 3D illustration
- Easy to use menu, similar to Smart Phone Apps
- Large buttons, easy to operate when wearing gloves
- Self-guidance to users to deal with specific situations
- Centered and angled down for easy reach and viewing
- Optional: 21 CFR Part 11 Compliance



Single Piece Wall

- Easy to reach service fixtures and electrical outlets on sidewalls
- Large radius corners for easy cleaning



User-friendly Work Tray

- Largest useable area in the market
- Recessed to contain spillage
- Sloped perimeter for easy cleaning
- Large, easy to clean tray handle



Raised Arm Rest

- Prevent grille blocking
- Comfortable working posture
- Durable stainless steel construction



Ergonomic Work Zone

- 10° angle to optimize user comfort, reduce glare, and maximize reach into the work area
- Brightly illuminated with >1200 lux (111 ft. cd)
- Industry-leading dimmable LED for optimum work comfort
- Airtight seal port for cable/tube exit protected by a negative pressure side wall





Esco Labculture® G4 Class II Type A2 Biosafety Cabinet
Available in 3 feet, 4 feet, 5 feet, and 6 feet models.

NET, FEATURING ADVANCED TOUCHSCREEN CONTROLLER

Energy-efficient DC ECM Blower

- The leading energy-efficient Class II Type A2 Biosafety Cabinet in the world with 70% energy savings compared to AC motor
- Stable airflow despite building voltage fluctuations and filter loading
- Standby mode to further reduce power consumption by 80%



Advanced ULPA Filtration System

- 10x Filtration efficiency of HEPA filter
- Creates ISO Class 3 work zone instead of industry-standard ISO Class 5
- Same 10 years filter life and replacement cost as HEPA filters




Dimmable LED Lamp

- Save energy and optimize work comfort




Removable Paper Catch

- Prevent objects from being pulled into blower plenum
- Removable for easy cleaning
- Optional pre-filter can be fitted



Tray Support Beams

- Support work tray evenly for less vibration
- Cleaning holder to easily wipe the drain pan



ISOCIDE™ Powder Coat

- Silver-ion impregnated powder coat
- Inhibits microbial growth to improve safety
- Prevents the plenum from becoming biohazard landfill


Certification


| | Performance | Air Quality | Filtration | Electrical Safety |
|----------------------|--------------------|--|---|---|
| Standards Compliance | NSF / ANSI 49, USA | ISO 14644.1, Class 3, Worldwide JIS B9920, Class 3, Japan JIS B55295, Class 3, Japan US Fed Std 209E, Class 1 USA | EN-1822 (H14), Europe IEST-RP-CC001.3, USA IEST-RP-CC007, USA IEST-RP-CC034.1, USA | UL 61010-1 3rd Ed, USA CSA22.2, No.1010-192, Canada EN-61010-1, Europe IEC61010-1, Worldwide |

Dynamic Chamber™

- Blower plenum and side walls are surrounded by negative pressure
- Prevent contaminants from escaping outside

- Positive Pressure
- Negative Pressure






Dynamic air barrier, where inflow and downflow converge
Side capture zones

- ULPA-filtered air
- Unfiltered / potentially contaminated air
- Room air / Inflow air

Cabinet Filtration System




- Ambient air is pulled through front grille to create inflow, without going into the work surface. Inflow is joined by half of the downflow, to create front air curtain that is fine-tuned to create a large performance envelope. The combined air stream travels through the back air column towards the blower.
- Approximately 1/3 of the air in the common plenum is exhausted through the ULPA filter to the room. The remaining 2/3 of the air is passed through the downflow ULPA filter and into the work area as a vertical laminar flow air to create ISO Class 3 work surface and prevents cross contamination.
- Near the work surface, the downflow splits. About half goes to the front grille, and half goes to the rear grille. A small portion enters the side capture zones to prevent dead air corners (small blue arrows).
- The design was optimized to give large performance envelope, that provides operator and product protection at wide Inflow and Downflow variation from the Nominal point.

The Performance Envelope Concept



- Nominal Airflow
- Personnel / Product Protection
- Area of no Personnel / Product Protection
- ▲ No Personnel / Product Protection
- Area of Personnel / Product Protection

Engineering Drawing



1. Exhaust sensor
2. USB Port
3. Zero Volt Relay Contact
4. Power Inlet
5. Centurion Touchscreen Controller

6. 10° Angled Sash Window
7. Max. Sash Opening
8. Single-piece Stainless Steel Back Wall
9. Single-piece Stainless Steel Work Tray
10. Electrical Panel

11. LED Lamp
12. IV Bar Retrofit Kit Provision
13. Service Fixture Retrofit Kit Provision
14. Duplex GFCI Outlet (L/R), Max 5 Amp total
15. Stainless Steel Arm Rest

16. Exhaust Collar (optional)
17. Exhaust Filter
18. DC ECM Blower
19. Downflow Filter
20. UV Lamp Provision
21. Cable Port

Options and Accessories

| | | | | | |
|--------------------------------------|------------------------------------|-------------------------------|-------------------------------|-------------------------------|-------------------------------|
| Anti-blowback Valve 10 inches | EG Powder Coated | ABBV-10P 5170352 | | | |
| | 304 Stainless Steel | ABBV-10S 5170354 | | | |
| Exhaust Collar | | ECO-F-LA2/AC2-3 G4 5171097 | ECO-F-LA2/AC2-4 G4 5171098 | ECO-F-LA2/AC2-5 G4 5171099 | ECO-F-LA2/AC2-6 G4 5171100 |
| UV Lamp | | UV-15A-L 5170251 | UV-30A-L 5170255 | | |
| IV Bar | | IV-955 5170276 | IV-1260 5170277 | IV-1565 5170278 | IV-1870 5170279 |
| Electrical Outlet | Direct Mounted | EO-H_ | | | |
| | GFCI | EO-GFCI 5170071 | | | |
| Service Fixtures | EU SF-Gas-20 mm and Solenoid Valve | SF-1G20 5170410 | | | |
| | EU SF-Vacuum-20 mm | SF-1V20 5170457 | | | |
| | EU SF-Air-20 mm | SF-1A20 5170502 | | | |
| | EU SF-Nitrogen-20 mm | SF-1N20 5170503 | | | |
| | EU SF-Water-20 mm | SF-1W20 5170458 | | | |
| | US SF-Universal-20 mm | SF-2U22 5170504 | | | |
| Support Stand | | STA-3 5131340 | STA-4 5131341 | STA-5 5131427 | STA-6 5131389 |
| Pipette Storage Shelf | | 5260327 | | | |
| Arm Rest Padding | | MEWREST 5170127 | | | |
| Foot Rest | | FT-REST 5170492 | | | |
| Laboratory Chair | | ME-LD-AR360 1150006 | | | |
| IQ OQ Protocol | | 9010179 | | | |




ABBV_



ECO-L_2_-MK3-LH



UV_A-L



IV_



EO-H_



EO-GFCI



SF-1_



SF-2U_



STA_



Pipette StorageShelf



MEWREST



FT-REST



ME-LD-AR360



IQQQ

TECHNICAL SPECIFICATIONS (10" OPENING)

| Labculture® Class II | Stainless Steel Side Walls | 110-130 VAC, 50/60 Hz | LA2-359 G4 10" 2011683 | LA2-459 G4 10" 2011685 | LA2-559 G4 10" 2011687 | LA2-659 G4 10" 2011689 |
|--|---|---|---|--|--|--|
| Nominal Size | | | 0.9 meter (3') | 1.2 meter (4') | 1.5 meter (5') | 1.8 meter (6') |
| External Dimensions* (W x D x H) | Without Arm Rest | | 1035 x 815 x 1570 mm (40.7" x 32.1" x 61.8") | 1340 x 815 x 1570 mm (52.8" x 32.1" x 61.8") | 1645 x 815 x 1570 mm (64.8" x 32.1" x 61.8") | 1950 x 815 x 1570 mm (76.8" x 32.1" x 61.8") |
| | With Arm Rest | | 1035 x 873 x 1570 mm (44.0" x 34.4" x 61.8") | 1340 x 873 x 1570 mm (52.8" x 34.4" x 61.8") | 1645 x 873 x 1570 mm (64.8" x 34.4" x 61.8") | 1950 x 873 x 1570 mm (76.8" x 34.4" x 61.8") |
| Internal Dimensions (W x D x H) | | | 915 x 625 x 720 mm (36.0" x 24.6" x 28.3") | 1220 x 625 x 720 mm (48.0" x 24.6" x 28.3") | 1525 x 625 x 720 mm (60.0" x 24.6" x 28.3") | 1830 x 625 x 720 mm (72.0" x 24.6" x 28.3") |
| Usable Work Area | | | 0.47 m ² (5.0 sq.ft.) | 0.62 m ² (6.7 sq.ft.) | 0.80 m ² (8.5 sq.ft.) | 0.96 m ² (10.3 sq.ft.) |
| Maximum Sash Opening | | | 580 mm (22.8") | | | |
| Average Airflow Velocity | Inflow | | 0.53 m/s (105 fpm) | 0.53 m/s (105 fpm) | 0.53 m/s (105 fpm) | 0.53 m/s (105 fpm) |
| | Downflow | | 0.30 m/s (60 fpm) | 0.30 m/s (60 fpm) | 0.33 m/s (65 fpm) | 0.30 m/s (60 fpm) |
| Airflow Volume | Inflow | | 446 m ³ /h (263 cfm) | 591 m ³ /h (350 cfm) | 741 m ³ /h (438 cfm) | 887 m ³ /h (525 cfm) |
| | Downflow | | 639 m ³ /h (374 cfm) | 848 m ³ /h (499 cfm) | 1063 m ³ /h (624 cfm) | 1272 m ³ /h (748 cfm) |
| | Exhaust | | 446 m ³ /h (263 cfm) | 591 m ³ /h (350 cfm) | 741 m ³ /h (438 cfm) | 887 m ³ /h (525 cfm) |
| | Required Exhaust with Optional Thimble Exhaust Collar | | 635 m ³ /h (374 cfm) | 627 m ³ /h (369 cfm) | 908 m ³ /h (534 cfm) | 1245 m ³ /h (733 cfm) |
| | Static Pressure for Optional Thimble Collar | | 33 Pa / 0.13 in H ₂ O | 18 Pa / 0.06 in H ₂ O | 37 Pa / 0.15 in H ₂ O | 61 pa / 0.25 in H ₂ O |
| ULPA Filter Typical Efficiency | | | ≥99.999% for particle size between 0.1 to 0.3 microns | | | |
| Sound Emission (dBA)* | NSF / ANSI 49 (Anechoic chamber) | | 56 | 59 | 60 | 61 |
| Light Intensity | | | > 1200 lux (111 ft. cd) | | | |
| Cabinet Construction | Main body | Electro-galvanized steel with white oven-baked epoxy-polyester Isocide™ antimicrobial powder-coated finish, 1.5 mm (0.06") / 16 gauge thick | | | | |
| | Work Zone | Stainless steel Type 304 with No.4 finish, 1.5 mm (0.06") / 16 gauge thick | | | | |
| Electrical | Nominal power | | 186 W | 200 W | 277 W | 385 W |
| | Heat Load | | 635 BTU/hr | 682 BTU/Hr | 945 BTU/Hr | 1314 BTU/Hr |
| | Full Load Amps 110-130V exclude 5A EO | | 10 A | | 15 A (2 power inlet) | |
| | Optional Outlets FLA | | 5A | | 5A | |
| Net Weight | | | 243 Kg (536 lbs) | 287 Kg (633 lbs) | 381 Kg (840 lbs) | 400 kg (882 lbs) |
| Shipping Weight | | | 292 Kg (644 lbs) | 350 Kg (772 lbs) | 439 Kg (968 lbs) | 506 kg (1116 lbs) |
| Shipping Dimensions, Maximum (W x D x H) | | | 1185 x 890 x 1900 mm (46.7" x 35.0" x 74.8") | 1490 x 890 x 1900 mm (58.7" x 35.0" x 74.8") | 1795 x 890 x 1900 mm (70.7" x 35.0" x 74.8") | 2100 x 890 x 1900 mm (82.7" x 35.0" x 74.8") |
| Shipping Volume Dimensions (W x D x H) | | | 2.00 m ³ (cu.ft.) | 2.5 m ³ (cu.ft.) | 3.00 m ³ (cu.ft.) | 3.6 m ³ (cu.ft.) |

*Noise reading in open field condition / anechoic chamber. Noise reading in normal room varies by room size, layout, and background noise, but may reach roughly 3-4 dBA above these values.

Life Science Tools



Biological Safety Cabinets



Laminar Flow Cabinet



Animal Workstations



CO₂/O₂ Incubators

Chemical Research



Ducted Fume Hoods



Ductless Fume Hoods



USP 797 / 800
Compliant Isolators



Powder Containment
Workstations

Pharma Compounding

Medical (ART)



Multi-room Incubators



Time-Lapse Incubators



IVF Workstations



PCR Workstations

PCR

Bioprocessing

Pharmaceutical



Bioreactors



cGMP Isolators



Filling Line Isolators



Downflow Booths

ESCO LIFESCIENCES GROUP

42 LOCATIONS IN 21 COUNTRIES ALL OVER THE WORLD



- 📍 Global Offices
- 📍 Joint Ventures
- 📍 Licensee
- 📍 Distributors
- 📍 Factories
- 📍 R&D Centers
- 📍 Regional Distribution Centers

Follow us on social media, download our apps,
and scan the QR code for more info.



@EscoLifesciences



@EscoLifesciences



@EscoLifesci



@Esco



@EscoLifesciences



@EscoLifesciences



Esco Lifesciences



Esco Lifesciences

ESCO[®]

LIFESCIENCES GROUP

Esco Technologies, Inc. • 903 Sheehy Drive, Suite F, Horsham, PA 19044, USA
Tel: +1 215-441-9661 • eti.admin@escolifesciences.com

Esco Micro Pte. Ltd. • 21 Changi South Street 1 • Singapore 486 777
Tel +65 6542 0833 • mail@escolifesciences.com
www.escolifesciences.com

Esco Lifesciences Group Offices: Bangladesh | China | Denmark | Germany | Hong Kong | India | Indonesia | Italy | Japan | Lithuania | Malaysia | Myanmar | Philippines | Russia | Singapore | South Africa | South Korea | Taiwan | Thailand | UAE | UK | USA | Vietnam